USER MANUAL

Hereby, MOB, declares that item MO2623 is in Compliance with the essential requirements and other relevant conditions of Regulation (EU) 2023/988. The full text of the EU declaration of conformity is available at the following internet address: www.momanual.com.

MOB, PO BOX 644, 6710 BP (NL).

PO: XXXXXXXX

Made in China

MO2623

EN

Adjustable relief cap

Instructions:

1. For cold compress:

Freeze the cap in a refrigerator for 2 hours. Then, take it out and put it on to relieve your pain or facial tension.

2. For hot compress:

Heat the cap in a microwave starting at 20 seconds and add in 10-second increments to get the desired temperature. Then, take it out and put it on to relieve your pain or facial tension. All microwaves may vary.

3. Do not wash in water. Wipe it with a damp towel and then allow to air dry.

Cautions:

- Before every use, thoroughly check the inner lining of the cap. Do not use if any gel is visible.
- Always freeze in a freezer bag. If it is too cold immediately out of the refrigerator, let it sit out until it reaches your desired temperature before placing comfortably on your head.
- Do not microwave for more than the time mentioned above. Do not apply if the cap is too hot to touch.
- Apply to dry skin only. Do not apply to infected areas or open wounds.
- If the gel comes in contact with any skin or eyes, immediately wash the area with soap and water. Consult your doctor if you still feel uncomfortable.

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